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N

Needle-Eze-extension	18
Needle-free self refilling syringe	12

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Syringes, needles and dosage equipment

A. Jeringas / Syringes

A-1. Jeringas manuales / Manual Syringes

Jeringa plástico 10 ml lisa

Plastic syringe, smooth non graduated 10 ml

Ref. 15.110

Jeringa plástico 20 ml lisa

Plastic syringe, smooth non graduated 20 ml

Ref. 15.112

Jeringa plástico 30 ml lisa

Plastic syringe, smooth non graduated 30 ml

Ref. 15.113

Jeringa plástico 50 ml lisa

Plastic syringe, smooth non graduated 50 ml

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Dermoject self-refilling syringe

Dedicated to intracutaneous injections for local anesthesia, corticoids injection and vaccinations

Regulated for a monodoses of 0,1 ml

The 4 ml of volume provides a capacity of 40 dose

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Fitted with a swiveling bottle support so that it can be used at any angle.

Typical applications in Veterinary Medicine are vaccination against Rabbit's Myxomatosis, vaccination meat pigs against Aujeszky's disease, induction of female Rabbits ovulation submitted to artificial insemination.

Produces a characteristic papule on the top of which can be seen a small hole where the liquid penetrated the integument.

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
SPARE PARTS FOR SOCOREX SELF-REFILLING SYRINGE

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
Needle-Free Injection System for Swine

14

The Pulse™ 250 is a highly flexible device designed to fit your current vaccination program. Designed to propel a jet of medication through the skin and into the muscle tissue, the new Pulse 250 is proven to increase profitability by eliminating the risks and expense associated with needles.

- *Eliminates needle fragments in meat*
- *Eliminates disease transmission from contaminated needles*
- *Sets new standard for worker safety*
- *Reduces injection site tissue damage*
- *Reduces vaccination costs*
- *Improves overall efficiency*

Handpiece

- *Handpiece manufactured from durable, polycarbonate medical grade material.*
- *Adapted to left- or right-handed operation, ergonomically designed to minimize hand fatigue.*
- *Actuated by contact pressure with pig, with no triggers to squeeze.*
- *Redundant on/off switch and hand-activated “enable” safety switch.*
- *High-pressure, Kevlar®-lined hose.*
- *Jeweled orifice assures optimal stream focus.*
- *Serrated nozzle head assures firm, “no slip” contact with animal during injection.*

Pneumatic Amplifier

- *Environmentally sealed Pneumatic Amplifier boosts, medication pressure at nozzle head.*
- *Offers multiple pressure settings to accommodate varying animal sizes.*
- *All components coming in contact with medication made of medical grade material for easy cleaning and maintenance.*

Energy System Vest

- *Fully adjustable, machine-washable Energy System | vest assures comfortable fit for all users.*
- *Chest pockets accommodate medicine bottles of any size, markers, gloves and other materials.*
- *Hip pockets contain CO2 bottle with regulator and amplifier.*
- *Each 24-ounce CO2 canister supplies power to inject 300 to 350 pigs.*
- *Coiled expansion hose.*
- *“Quick disconnect” attachments for easy removal, cleaning adjustment.*

Directions for Use

1. Preparing the system for use.

Pulse 250 comes complete with a User Passport. Begin each injection sequence by recording the beginning injection number shown on the counter mounted in the Pneumatic Amplifier, the type of medication administered, animal identity numbers or other data. It records ending count at conclusion of each injection sequence.

2. Fitting and adjusting the Energy System.

The fully adjustable vest is designed to assure comfortable fit and weight distribution for users of all sizes and statures. Simply connect and adjust the two quick-snap chest straps and waist belt, insert the Pneumatic Amplifier and fresh CO2 cartridge in the side pockets, and adjust system so it is comfortably positioned.

3. Pressuring the system.

Firmly connect the threaded end of the regulator to the top of the CO2 cartridge, and attach the quick-connect coupler to the Pneumatic Amplifier. Insert the needle of a standard draw-off into the medicine bottle, and attach draw-off to the medicine chamber on the Pneumatic Amplifier. Turn CO2 valve to the full-open position.

4. Select proper pressure setting.

A regulator located on the Pneumatic Amplifier offers multiple pressure settings to accommodate the size of animal you will be injecting. After setting the dial to the appropriate pressure setting, test the system by administering four or five test injections into the foam test pad provided with your Pulse 250.

5. Begin injection operations.

Designed to recycle approximately every one second, Pulse 250 allows you to inject the smallest nursery pigs as rapidly as conditions permit.

6. Press nozzle firmly against the skin.

Pulse 250 is operated through firm nozzle contact with the animal. Place handpiece nozzle firmly against the desired injection site. Press lightly (about 1-1/2 pounds offered), and the pressurized system will automatically inject a 1 to 2,5 milliliter dose of medicine.

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Dosatron D200RE

Characteristics:

Dosage technique: Proportional without electricity

Energy source: The pressure of the water

Dosage values: p.e. 1%=1:100=1V concentrated product

Regulation of the dosage: External, modifiable during the operation

Value min. - max. of the flow of water: 20-1505 l/h

Operation pressure: 0.5 to 5 Bar (entrance pressure)

Cylinder capacity: 0.42 l (1 cycle)

Operation temperature: max. 40°C, min. protect from freezes

Viscosity max. of the concentrate: 400 cst (at 20°C). Kit PV recommended starting from 200 cst

Ref. 24.101

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Mixer agitation set A125 12V.

Ref. 24.101.01